



- Summary Report -


Test Product:	Tesa + Ipa + CJC 6/3/3mg	Manufacturer:	N/A
Client:	Lumineris Scientific	Batch/Lot #:	LS-TiC12-2681
Description:	3ml glass vial of blue, lyophilized powder. (Trans Green Cap / Purple crimp)	Receive Date:	08-13-26

COMPONENT	TEST	SPECIFICATION	RESULT
1. Tesamorelin	Quantitative Assay	6mg/vial (NLT 90%)	6.45mg ∴ Pass
2. Ipamorelin	Quantitative Assay	3mg/vial (NLT 90%)	3.46mg ∴ Pass
3. CJC-1295	Quantitative Assay	3mg/vial (NLT 90%)	3.60mg ∴ Pass
	Percent Purity	NLT 95%	99%
	Correlation Coefficient	$r_{xy} = \frac{\sum_{i=1}^n (x_i - \bar{x}) \cdot (y_i - \bar{y})}{\sqrt{\sum_{i=1}^n (x_i - \bar{x})^2} \cdot \sqrt{\sum_{i=1}^n (y_i - \bar{y})^2}}$	∴ Pass

Tested on a Jenway 6715
 UV/Vis Spectrophotometer.

Discussion:

Sample received sealed
 and in good condition.



Performing Analyst Initial & Date:  23 August 2026

This data is presented in good faith and not warranted as to the accuracy of results. Use of this information "as is" shall be at the risk of the user, and independent verification is the sole responsibility of the user. Analytical Formulations Inc. disclaims any liability incurred related to the misuse of this data. This report does not grant permission to assert any nutritional claims, nor does it give license to use any existing patents or copyrights pertaining to the product samples, or any individual component of said goods. This report may not be reproduced, except in its totality, without the prior, express written consent of Analytical Formulations Inc.